

Thin Film Stress Measurement at NUFAB

Ying Jia

NUFAB
NORTHWESTERN UNIVERSITY
MICRO/NANO FABRICATION FACILITY

NUANCE
Northwestern University Atomic and
Nanoscale Characterization Experimental Center

Northwestern
EXPLORING INNER SPACE

Quick Survey

Did you ever consider film stress in film deposition process or device fabrication or film characterization?


Yes/No in chat

Outline

- Where does the film stress come from?
- Why should you care about the stress?
- How to measure the stress?
- The Toho-FLX 2320-S thin film stress measurement system
- Control of the stress


Where does the stress come from?

- Induced during film deposition
- Intrinsic stress:
Non-equilibrium nature of deposition
Lattice mismatch, impurities etc.
- Extrinsic stress:
Environment change
Thermal expansion coefficients mismatch
Nonuniform plastic deformation




Thermal stresses

- In a structure with inhomogeneous thermal expansion coefficients subjected to a uniform temperature variation
- In a homogeneous material exposed to a thermal gradient





Consequences of high stresses

The adhesion is not very good




The adhesion is very good




Stress measuring techniques

The deflection technique

- A stressed thin film will bend a moderate thick substrate by a measurable degree
- Measure the curvature or deflection of the substrate before and after coating
- Simple and fast



Calculate stress



$$\delta(x) = a + bx + cx^2$$

Radius


$$R = \frac{1}{2c}$$

$$\sigma_f = \frac{E_s}{(1 - \nu_s)} \frac{t_s^2}{t_f} \left(\frac{1}{R_a} - \frac{1}{R_b} \right)$$

E_s, ν_s : Young's modulus, Poisson ratio of substrate
 t_s, t_f : thickness of substrate and film
 R_a, R_b : radius before and after coating

Film stress measurement at NUFAB

Toho FLX 2320-S




- Two laser (670 nm and 785 nm) to resolve the possible destructive interference
- Measurement Range 1 MPa to 4 GPa
- Accuracy Less than 2.5% or 1 MPa (whichever is larger)
- Scan range programmable up to 200mm
- Minimum scan step 0.02 mm
- 3D mapping
- In-situ stress measurements from room temperature to 500°C

Limitations:


1. Not local stress
2. Too much roughness, low reflectivity
3. Transparent substrate – use Dektak

3D mapping



Control of the stress

- Adjust coating process parameters
 - Film thickness
 - Deposition temperatures
 - Deposition rate
 - Substrate etc..





J. Vac. Sci. Technol. A **36**, 020801 (2018)

- Thermal relaxation
 - Crystalline slip

Toho FLX 2320-S

In-situ stress measurements from room temperature to 500°C





Contact me: ying.jia@northwestern.edu

General NUFAB contact: nufab@northwestern.edu